



HUBER NEEDLE RANGE

POLYPERF® - Huber needle with connecting line

HUBER NEEDLES

- ▶ Complete range of sizes and gauges
- ▶ For the puncture of venous, intra-rachidian, arterial and intra-peritoneal implantable ports

References	Description	Usable length (mm)	Recommended field of use	Packaging
522506	Curved	6/10 ^e 23G	20	Box of 50
522207		7/10 ^e 22G	25	
523807		7/10 ^e 22G	35	
522509		9/10 ^e 20G	20	
523809		9/10 ^e 20G	35	
522511		11/10 ^e 19G	20	
523811		11/10 ^e 19G	35	
512506	Straight	6/10 ^e 23G	25	Box of 50
512507		7/10 ^e 22G	30	
513807		7/10 ^e 22G	38	
512509		9/10 ^e 20G	25	
512511		11/10 ^e 19G	25	
513811		11/10 ^e 19G	38	



POLYPERF®

- ▶ Huber needle with connecting line - **DEHP FREE (without phthalates) – LATEX FREE**
- ▶ No interaction with anti-cancerous agents

POLYPERF® without lateral injection site – Huber needle

References	Description	Usable length (mm)	Tubular (cm)	Packaging
581707	Curved needle with luer lock and clamp	7/10 ^e mm 22G	17	Box of 10
582007		7/10 ^e mm 22G	20	
582507		7/10 ^e mm 22G	25	
583007		7/10 ^e mm 22G	30	
583507		7/10 ^e mm 22G	35	
581709		9/10 ^e mm 20G	17	
582009		9/10 ^e mm 20G	20	
582509		9/10 ^e mm 20G	25	
583009		9/10 ^e mm 20G	30	
583509		9/10 ^e mm 20G	35	
581711		11/10 ^e mm 19G	17	
582011		11/10 ^e mm 19G	20	
582511		11/10 ^e mm 19G	25	
583011		11/10 ^e mm 19G	30	
583511		11/10 ^e mm 19G	35	



POLYPERF® with lateral injection site – Huber needle

References	Description	Usable length (mm)	Tubular (cm)	Packaging
592007	Curved needle with luer lock and clamp with lateral injection site	7/10 ^e mm 22G	20	Box of 10
592507		7/10 ^e mm 22G	25	
592009		9/10 ^e mm 20G	20	
592509		9/10 ^e mm 20G	25	
592011		11/10 ^e mm 19G	20	
592511		11/10 ^e mm 19G	25	

